SPECTROLINE® SB-100P Series Long Wave UV Lamps

Ultraviolet radiation emitted from this product. Avoid exposure. ALWAYS WEAR PROTECTIVE CLOTHING. EXPOSURE MAY CAUSE PREMATURE AGING OF THE SKIN AND CANCER. ALWAYS WEAR PROTECTIVE EYEWEAR; FAILURE TO DO SO MAY RESULT IN SEVERE BURNS OR LONG TERM INJURY TO EYE. Never look directly into the lamp. Exposure can cause eye and skin allergy and allergic reactions. Medications or cosmetics may increase your sensitivity to ultraviolet radiation. Consult physician before operating this product if you are using medications or have a history of skin problems or believe yourself especially sensitive to sunlight.

FOR PROFESSIONAL EXAMINATION USE ONLY

IMPORTANT SAFETY INSTRUCTIONS

1. Read all instructions.
2. Care must be taken, as burns can occur from touching hot parts.
3. Do not operate the lamp with a damaged cord or if the lamp has been damaged.
4. **Do not look directly into the UV light. Wear UV-absorbing protective eyewear (such as the UVS-30 glasses provided) when operating this or any UV lamp. Never operate the lamp with the UV filter lens removed.**
5. Never operate this equipment in any manner not specified in these instructions because the protection provided by the equipment may be impaired.
6. Adequate ventilation should be provided when working on operating combustion engines.
7. This unit is not recommended for high-dust environments.

GENERAL INFORMATION

1. Install the UV bulb and the Bulb-Saver™ (packed separately).
2. Connect the plug to a power outlet that complies with the electrical requirements specified on the lamp's label.
3. Turn ON the lamp. Allow approximately 3 to 5 minutes for the lamp to reach full lighting intensity. If the lamp fails to light, check that the bulb is completely screwed into the socket.
4. When the lamp is used intermittently and frequently, do not turn it off until the inspection is complete. Any interruption in the power supply, such as a drop in voltage or changing power outlets, will cause the bulb to go out. The bulb must then cool for a minimum of 5 minutes before restarting. Continuous operation will also extend the bulb life.
5. The lamp housing and filter will become hot during continuous operation. Do not touch these components. We recommend using the *optional* Spectroline HGS-100A heat guard/stand for increased protection.

OPERATOR SERVICING

Do not perform any servicing that is not covered in these instructions. Return the lamp to the factory for maintenance not covered in this section.

HOUSING AND FILTER MAINTENANCE

- **Always unplug the lamp from the power outlet before performing any maintenance or servicing. Allow the UV lamp to cool before handling.**
- Immediately clean spills from the housing with water and a mild detergent and wipe dry. Periodically clean the filter with a glass cleaner and a soft cloth.

BULB/FILTER INSTALLATION OR REPLACEMENT

Use only a Spectroline® 100S/M UV bulb. Using any other bulb may damage the unit and will void the warranty.

1. Detach the silicone-rubber filter holder by peeling it from around the circumference of the lamp housing. The filter will remain inside the filter holder. (If replacing the filter, press it firmly but gently inside the filter holder.) Set the filter holder aside.
2. Place the UV bulb through the Bulb-Saver™/collar in the lamp housing. Turn the bulb *clockwise* gently but firmly into the socket with the fingertips. (If removing an existing bulb, turn it *counterclockwise*.) **Do not operate the lamp without the Bulb-Saver.**
3. Reattach the filter holder by pressing it with the heel of your hand around the circumference of the lamp housing. Be sure the groove in the filter holder seats over the rim of the housing.

UV INTENSITY VERIFICATION

To ensure that your SB-100P is operating at the required UV intensity, it should be checked periodically. The AccuMAX™ XRP-3000 radiometer/photometer kit uses a dual-wavelength UV-A/VIS sensor detector to accurately measure ultraviolet irradiance or visible light. The AccuMAX XR-1000 digital readout unit with the XS-365 UV sensor or the Spectroline DM-365XA digital readout unit are also recommended for accurate UV-A measurement. The meters are specially designed to measure the UV irradiance from 320–400nm with a peak at 365nm.

For more information about these units, please contact the Customer Service Department at Spectronics Corporation. In the U.S. and Canada, call toll-free 1-800-274-8888. Outside the U.S. and Canada, call 516-333-4840.

SPECIAL PURPOSE FILTERS

- The SB-100PC Series lamps are equipped with a 2F100C clear filter that reduces transmission of medium wave (UV-B) and short-wave (UV-C) radiation to safe levels, while permitting high UV-A energy for curing applications.
- The SP-100PY semiconductor wafer inspection lamp is equipped with a 2F475Y yellow filter that blocks the emission of wavelengths shorter than 500nm, protecting against activation of the photoresist during semiconductor wafer inspection.
- The SB-100PD Series lamps are equipped with a 2F350 diffusing filter that produces a "diffused spot" beam of UV-A light.

ELECTRICAL SPECIFICATIONS

Model	Filter	120V/60Hz	230V/50Hz	240V/50Hz	100V/50-60Hz
SB-100P	2F958	1.05 Amps	—	—	—
SB-100PD	2F350	1.05 Amps	—	—	—
SB-100PC	2F100C	1.05 Amps	—	—	—
SB-100PY	2F475Y	1.05 Amps	—	—	—
SB-100P/F	2F958	—	0.95 Amps	—	—
SB-100PD/F	2F350	—	0.95 Amps	—	—
SB-100PC/F	2F100C	—	0.95 Amps	—	—
SB-100PY/F	2F475Y	—	0.95 Amps	—	—
SB-100P/FA	2F958	—	—	0.95 Amps	—
SB-100PD/FA	2F350	—	—	0.95 Amps	—
SB-100PC/FA	2F100C	—	—	0.95 Amps	—
SB-100PY/FA	2F475Y	—	—	0.95 Amps	—
SB-100P/FB	2F958	—	0.95 Amps	—	—
SB-100PD/FB	2F350	—	0.95 Amps	—	—
SB-100PC/FB	2F100C	—	0.95 Amps	—	—
SB-100PY/FB	2F475Y	—	0.95 Amps	—	—
SB-100P/J	2F958	—	—	—	2.7 Amps
SB-100PD/J	2F350	—	—	—	2.7 Amps
SB-100PC/J	2F100C	—	—	—	2.7 Amps
SB-100PY/J	2F475Y	—	—	—	2.7 Amps

All bulbs are 100S/M, 100W Medium Base Spot

ENVIRONMENTAL CONDITIONS

The SB-100P Series UV lamps are designed to be safe under the following conditions:

- Indoor use;
- Altitude up to 2,000 m (6,562 ft.);
- Temperature 5°C to 40°C (41°F to 104°F);
- Maximum relative humidity 80% for temperatures up to 31°C (88°F) decreasing linearly to 50% relative humidity at 40°C (104°F);
- Mains supply voltage fluctuations not to exceed ±10% of the nominal voltage;
- Installation Category II;
- Pollution Degree 2.

WARRANTY

The warranty policy for the SB-100P Series UV lamps is provided on the Certificate of Limited Warranty enclosed separately with each unit.

NOTE: For assistance of any kind, contact the Customer Service Department at Spectronics Corporation.

Item	Part #	Description
1	123604	Filter Holder
2	2F958	Filter Lens, UV
	2F350	Filter Lens, UV, Diffusing
	2F100C	Filter Lens, Clear
	2F475Y	Filter Lens, Yellow
3	100S/M	Bulb, 100W, Medium Base Spot
4	106862	Bulb-Saver™
5	EC-0116	Socket, Medium Base, 120V
	102931	Socket, Medium Base, 230/240V
	EC-0116	Socket, Medium Base, 100V
6	124736	Lamp Housing
7	103767	Plate, Rear Cover
8	100310	Pistol Grip
9	118677	Handle, Ballast
10	MC-0157	Mounting Block
11	MC-0158	Pin
12	SA-0101	Ballast Assembly, 120V
	105049	Ballast Assembly, 230/240V
	118989	Ballast Assembly, 100V
13	119799	Cord, Secondary, 8-ft
	119800	Cord, Secondary, 20-ft
14	119797	Cord, Primary, 8-ft*
	119798	Cord, Primary, 20-ft*
15	EC-0358	Power Switch

*Each of the primary cords is fitted with the appropriate plug for the country of shipment: Domestic, 120V; European, 230V; British, 230V; Australian, 240V or Japanese, 100V. For additional information, contact Spectronics Corporation.

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HINWEIS: Wenden Sie sich zwecks etwaiger Hilfestellung bitte an die Kundendienstabteilung der Spectronics Corporation.